

DIYNAFLUOR: Protocol for Low-Cost 0-10 ng/μL DNA Quantification Assay

Version 0.1

6/12/2024

Introduction

This low-cost DNA quantification assay costs ~\$0.13 USD per measurement and has a working range of ~0.3 – 12.5 ng/μL, with an average absolute bias (accuracy) of 0.218 ± 0.149 ng/μL.

For more information, please refer to our Preprint “*DIYNAFLUOR: An Affordable DIY Plug-and-Play Nucleic Acid Fluorometer for eDNA Quantification in Resource Limited Settings*”, which will be made available soon.

Materials

- *SYBR Safe 10,000X Concentrate* (Thermo Fisher Scientific S33102)
- *Lambda DNA 300 ng/μL* (Thermo Fisher Scientific SD0011)
- *LTE Buffer* (10 mM TRIS, 0.1 mM EDTA, pH 8.0) or 1 x TE Buffer (Thermo Fisher Scientific 12090015)
- *1.5 mL Microcentrifuge Tubes* (For preparing reagents and samples)
- *0.5 mL Thin-walled PCR Tubes* (We recommend Axygen PCR-05-C)
- *P1000, P200, P20 Disposable Pipette Tips*
- *Aluminium Foil*

Equipment

Required

- *P1000, P200, P20 Adjustable Pipettes*
- *Fine Point Pen* (For labelling Tubes)

Desirable

- *Tube Racks*
- *Vortex*
- *Mini Centrifuge*

Protocol

1. Prepare a 10 ng/μL DNA Standard.

We recommend preparing the 10 ng/μL DNA standard from 300 ng/μL Lambda DNA (Thermo Fisher Scientific), however any dsDNA with a known concentration should work.

For example, for 300 μL of 10 ng/μL Lambda DNA, add 10 μL of 300 ng/μL Lambda DNA to 290 μL of LTE Buffer.

Vortex for 2-3 seconds to mix and briefly centrifuge. (Alternatively, gently pipette mix 10-20x with the P200 set to 200 μL. Bring the contents to the bottom of the tube by gently tapping the bottom of the tube on a table).

Note: The 10 ng/μL DNA standard can be stored at -20 °C for future use, but repetitive freeze thawing should be avoided.

2. Prepare 200X SYBR Safe Stock Solution by diluting SYBR Safe 10,000X Concentrate in LTE.

For example, for 500 μL of 200X SYBR Safe Stock Solution, add 10 μL of SYBR Safe 10,000X Concentrate to 490 μL of LTE Buffer.

Vortex for 2-3 seconds to mix and briefly centrifuge. (Alternatively, gently pipette mix 10-20x with the P200 set to 200 μL. Bring the contents to the bottom of the tube by gently tapping the bottom of the tube on a table).

Keep away from light by wrapping the tube in Aluminium Foil.

The 10 ng/μL DNA Standard can be stored at -20 °C for future use (up to 4 months), but repetitive freeze thawing should be avoided.

3. Prepare SYBR Safe Master-Mix by adding 1 μL of 200X SYBR-Safe Stock Solution to 199 μL LTE for each sample and the two standards.

For example, to measure 5 Samples and 2 Standards, Add 7 x 1 μL = 7 μL of 200X SYBR-Safe Stock Solution to 7 x 199 μL = 1393 μL of LTE buffer.

4. Set up two 0.5 ml Thin-wall PCR Tubes for the Standards and one Tube for each Sample. Label the lids of the Tubes.

Label the Standard Tubes '1' and '2'. Label the Sample Tubes with names that allows them to be easily identified. Do not write labels on the sides of the Tubes.

5. Prepare the Standards.

Add 190 μL of the *SYBR Safe master-mix* to both *Standard Tubes 1 and 2*. Add 10 μL of *LTE* to *Standard Tube 1*. Add 10 μL of the *10 ng/ μL DNA Standard* to *Standard Tube 2*.

Vortex for 2-3 seconds to mix and briefly centrifuge. (Alternatively, gently pipette mix 10-20x with the P200 set to 200 μL . Bring the contents to the bottom of the tube by gently tapping the bottom of the tube on a table).

6. Prepare the Samples.

Add 190 μL of the *SYBR Safe master-mix* to all *Sample Tubes*. Add 10 μL of sample to all *Sample Tubes*.


Vortex for 2-3 seconds to mix and briefly centrifuge. (Alternatively, gently pipette mix 10-20x with the P200 set to 200 μL . Bring the contents to the bottom of the tube by gently tapping the bottom of the tube on a table).

7. Incubate all Tubes for 2 minutes at room temperature out of direct light.

8. Measure *Standard Tube 1*.

Insert *Standard Tube 1* into the DIYNAFLUOR. Ensure the “*High Sensitivity*” mode is selected, and the “*Sample Input*” is set to “10 μ L”.

Press “*Measure*”.




A data point for *Standard Tube 1*'s Fluorescence will be plotted.

Note: Ensure bubbles are not present in Tubes prior to measurements. Spin or gently tap the bottom of the tube on a table to remove them.

9. Measure *Standard Tube 2*.

Insert *Standard Tube 2* into the DIYNAFLUOR.

Press “*Measure*”.




A data point for *Standard Tube 2*'s Fluorescence and a 2-point calibration curve will be plotted.

10. Measure Samples.

Insert a *Sample Tubes* into the DIYNAFLUOR.

Press “*Measure*”.



A data point for the measured “*Tube Concentration*” will be plotted on the calibration curve. The “*Sample Concentration*” is the reported in lower right corner of the DIYNAFLUOR GUI. (The “*Sample Concentration*” is calculated from the dilution factor of Sample added to the *SYBR Safe master-mix*).


Repeat this process for the remaining samples.

Note: If a Sample’s measured Fluorescence is ~25% higher than Standard Tube 2’s Fluorescence, the Sample should be diluted and remeasured to fall within the 2-point calibration curve. If a Samples measured Fluorescence is within $\sim \pm 120$ units of Standard Tube 1’s Fluorescence value, then the Sample has insufficient DNA for measurement.

11. Save the Data.

Press “*Save*” to export the data as a .csv file.

Example Data: A 2-fold dilutions series of Lambda DNA measured from 0-10 ng/μL.



| | A | B | C | D | E |
|---|-------------|------------------------------|-----------------------------|----------------------------|-------------------|
| 1 | Name | Sample Concentration (ng/uL) | Measured Flouresence (arb.) | Tube Concentration (ng/uL) | Sample Input (uL) |
| 2 | Standard #0 | - | 1350.33 | 0 | - |
| 3 | Standard #1 | - | 4966.33 | 0.5 | - |
| 4 | Sample #0 | -0.224 | 1269.33 | -0.0112 | 10 |
| 5 | Sample #1 | 1.0998 | 1748 | 0.055 | 10 |
| 6 | Sample #2 | 2.5092 | 2257.67 | 0.1255 | 10 |
| 7 | Sample #3 | 5.106 | 3196.67 | 0.2553 | 10 |
| 8 | Sample #4 | 10.2913 | 5071.67 | 0.5146 | 10 |